



SOTA PT100T pH electrode, 12 mm, TOP68 connector

Proizvod #: LZX583
EUR Cijena: Kontaktirajte nas
Nije prikazan datum isporuke

Standard glass and reference system with head connector TOP68

Specifikacije

Dijafragma:	Ring diaphragm, porous PTFE, double
Dubina uranjanja:	105 mm
Duljina:	140 mm
Elektrolit:	KCl/AgCl + KNO ₃ gel, suitable for high temperatures
Garancija:	6 mjeseci
Konektori:	TOP68
Materijal:	Sensor: Glass
Mjerni rang:	0 - 80 °C
Navoj senzora:	Pg 13,5
Primjena:	General use
Rang tlaka:	0 - 3,5 bar
Referentna elektroda:	None
Spajanje kabla:	TOP68 connector
Temperaturni senzor:	PT100