



Bottle Carrier, Linear

Proizvod #:	2094300
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Linear, low-density polyethylene. Securely holds up to six standard 0.5 L glass bottles, pint-size reagent bottles, and most other 500 mL glass or plastic bottles—making breakage almost impossible. Lightweight, chemical-resistant, and easy to clean. Dimensions: 34.9 x 24.8 x 27.6 cm (14 x 10 x 11"). Maximum bottle size it will hold: 73 mm in diameter, 178 mm in height.

Specifikacije

Dubina:	11 Inches
Duljina:	349 mm
Širina:	10 Inches
Volumen/ veličina paketa:	1 each