



Be Right™



Scoop, sampler, sterile

Proizvod #:

2600900

EUR Cijena:

Kontaktirajte nas

Nije prikazan datum isporuke

Sterile scoop samplers hold various capacities. One-piece molded scoops used to obtain samples with one scoop. White, high-impact styrene. Individually bagged, lengths measuring from 17.8 to 21.6 cm (7 to 8.5"); 100 to case. Capacity: 60 mL (20 oz). Bowl: 25 x 76 x 38 mm (1 x 3 x 1.5") Bel-Art.

Specifikacije

Kapacitet: 60 mL

Volument/ veličina paketa: 1