



Thermometer, Non-Mercury, 0 - 220 °F

Proizvod #: 2676200
EUR Cijena: Kontaktirajte nas
Nije prikazan datum isporuke

Compact and durable, non-mercury thermometers enclosed in a protective nylon case. These thermometers are designed for ruggedness and come in many standard ranges. Suitable for all field and portable uses where an unprotected thermometer may break.

Specifikacije

Podjela termometra:	2 °F
Temperaturni rang:	0 °F - 220 °F
Tip termometra:	Non-Mercury
Volumen/ veličina paketa:	1 each