



Flask, volumetric, polypropylene, 100 mL

Proizvod #:	1406042
EUR Cijena:	Kontaktirajte nas
Nije prikazan datum isporuke	

Polypropylene flask and screw cap. Durable with excellent chemical resistance, with the same accuracy as glass. Permanent color graduation line with marking spot. Calibrated “to contain/deliver” at 20°C. The 200/203 capacity flask has dual marks useful in DO determination via Winkler Titration.

Specifikacije

Certificiranja:	Not Applicable
Jedinica:	1 each
Kapacitet:	100 mL
Materijal tijela mjerača:	Polypropylene