



Formaldehyde cuvette test - ISO 12460, 0.5-10 mg/L H₂CO, 25 tests

Proizvod #:

LCK425

EUR Cijena:

Kontaktirajte nas

Dostupno

LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer

Differentiation between tests and measuring ranges by means of colour-coding

Specifikacije

Broj testova:	25
Čuvanje uređaja:	15 - 25 °C
Mjerni rang:	0,5 - 10 mg/L H ₂ CO
Naziv metode:	ISO12460
Parametar:	Formaldehyde
Platforma:	LCK
Uređaj:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000